


Top DNA Polymerase

Top DNA polymerase is a novel thermostable polymerase that is more processive than *Taq*. In fact, the extension rate of Top is > 3 x that of *Taq* polymerase! Top polymerase can be used for a variety of PCR applications and is a robust enzyme for everyday PCR.



Applications

- Standard PCR up to 10 kb
- TA Cloning
- DNA Sequencing

Fig 1. Long PCR amplification test of Top DNA polymerase and *Taq* DNA polymerase using Lambda DNA. Note robust amplification from Top polymerase out to 10 kb.

Lane M1: 1 kb DNA Ladder (D-1040)
 Lane M2: Lambda DNA/Hind III Marker (D-1050)
 Lane 1: 2 kb PCR product
 Lane 2: 3 kb PCR product
 Lane 3: 4 kb PCR product
 Lane 4: 5 kb PCR product

Lane 5: 6 kb PCR product
 Lane 6: 7 kb PCR product
 Lane 7: 8 kb PCR product
 Lane 8: 9 kb PCR product
 Lane 9: 10 kb PCR product
 *10 ng of Lambda DNA and 1U of each DNA polymerase used for amplification.

HotStart DNA Polymerase

Bioneer's HotStart DNA polymerase uses an exclusive enzyme-mediated HotStart method that, unlike most other HotStart chemistries, completely releases all polymerase activity during the first denaturation step. Top DNA polymerase is completely inhibited by pyrophosphate at temperatures below 70°C. However, at temperatures above 70°C, a thermostable pyrophosphatase initiates pyrophosphate hydrolysis and activates the DNA polymerase. This prevents the formation of non-specific products and primer-dimers during the reaction set-up process and results in improved PCR specificity.

Applications

- HotStart PCR up to 12 kb
- PCR with complex genomic templates/low copy templates/cDNA
- Multiplex PCR reactions




Fig 2. Multiplex PCR comparison of genomic DNA using 6 sets of primers and 2 different DNA polymerases. Note uniform amplification of multiplex products with Bioneer's HotStart chemistry.

Lane M: 100 bp DNA Ladder (D-1030)
 Lane 1: 750 bp fragment
 Lane 2: 590 bp fragment
 Lane 3: 450 bp fragment


Lane 4: 360 bp fragment
 Lane 5: 260 bp fragment
 Lane 6: 150 bp fragment
 Lane 7: Multiplex PCR with primers used for Lane 1 to 6

TLA DNA Polymerase

AccuPower® TLA DNA polymerase provides both high fidelity and long range PCR capabilities. TLA DNA polymerase also features an enhanced processivity (TLA's base incorporation rate is almost 2 x as fast as *Taq*) that makes it ideal for your long range and standard PCR applications.

Applications

- Gene synthesis
- PCR or primer extension requiring high fidelity (up to 20 kb)
- Blunt end PCR cloning
- Mutagenesis



M: 1 kb DNA Ladder (D-1040)
 TLA: TLA DNA polymerase
Taq: *Taq* DNA polymerase
 Lane 1: 10 pg of λ DNA
 Lane 2: 1 pg of λ DNA
 Lane 3: 100 fg of λ DNA
 Lane 4: 10 fg of λ DNA

Fig 3. 1 kb amplification of Lambda DNA using TLA DNA polymerase and other competitors' proofreading enzymes. Note that TLA efficiently amplifies very low copy numbers or template.

CycleScript™ Reverse Transcriptase

CycleScript Reverse Transcriptase synthesizes cDNA from RNA with high activity across a wide range of temperatures from 37°C to 55°C. Reverse transcription can thus be carried out similar to PCR (Cyclic RT). The Cyclic RT reaction is composed of the following steps and is repeated to generate homogeneous cDNA: incubation at 15 to 40°C for primer annealing and extension; heating up to 45 to 48°C for extension (optional); and finally incubation at 50 to 55°C for denaturation of the secondary structure of the RNA. This kit can also be used for conventional reverse transcription.


Fig 4. Comparison of transferrin receptor gene amplification with Cyclic RT, conventional RT and other companies' products. Note specific and robust amplification of this rare gene using a Cyclic RT reaction with CycleScript.

Lane MW: 100 bp plus DNA Ladder (D-1035)
 Lane 1 - 4: TFR (Transferrin receptor gene) amplified with MMLV (E-3121)
 Lane 5 - 8: TFR amplified with CycleScript (E-3131)
 Lane 9 - 12: TFR amplified with CycleScript (E-3131)
 Lane 13 - 16: TFR amplified with MMLV from company I
 Lane 17 - 20: TFR amplified with S-script from company S
 Lane 21 - 24: TFR amplified with S-script-II from company I
 Lane 25 - 28: TFR amplified with S-script-III from company I
 Lane 29 - 32: TFR amplified with O-script from company Q

*700 ng of total RNA used for RT and the same amount of RT products used for electrophoresis.

Applications

- Gene synthesis
- First-strand cDNA synthesis from RNA, up to 9 kb
- cDNA or RNA sequencing
- Random primer probe labeling
- cDNA library construction
- RT-PCR
- mRNA 5'-end mapping by primer extension



© 2010 Bioneer Corporation. Bioneer® and Accubower® are registered trademarks of Bioneer Corporation. Prices and specifications are subject to change without notice. Order No. B10007.

Try before you buy!

Visit <http://us.bioneer.com/sample> and qualify for a free sample of our enzymes. In addition, if you respond to a short survey after you try our product, you can take 50% off your first order!

Ordering Information

Product Description	Cat No.	Price
Top DNA Polymerase, 500 U, 10 mM dNTPs, 10 x Reaction Buffer with MgCl ₂	E-3100	\$119.00
Top DNA Polymerase, 500 U, 10 mM dNTPs, 10 x Reaction Buffer, 20 mM MgCl ₂	E-3100-1	119.00
Top DNA Polymerase, 500 U, 10 x Reaction Buffer with MgCl ₂	E-3100-2	83.00
Top DNA Polymerase, 500 U, 10 x Reaction Buffer, 20 mM MgCl ₂	E-3100-3	83.00
Top DNA Polymerase, 2,000 U, 10 mM dNTPs, 10 x Reaction Buffer with MgCl ₂	E-3101	407.00
Top DNA Polymerase, 2,000 U, 10 mM dNTPs, 10 x Reaction Buffer, 20 mM MgCl ₂	E-3101-1	407.00
Top DNA Polymerase, 2,000 U, 10 x Reaction Buffer with MgCl ₂	E-3101-2	303.00
Top DNA Polymerase, 2,000 U, 10 x Reaction Buffer, 20 mM MgCl ₂	E-3101-3	303.00
HotStart DNA Polymerase, 250 U, 10 mM dNTPs, 10 x Reaction Buffer with MgCl ₂	E-3150	120.00
HotStart DNA Polymerase, 250 U, 10 x Reaction Buffer with MgCl ₂	E-3150-1	87.00
HotStart DNA Polymerase, 1,000 U, 10 mM dNTPs, 10 x Reaction Buffer with MgCl ₂	E-3151	384.00
TLA DNA Polymerase, 250 U, 10 mM dNTPs, 10 x Reaction Buffer with MgCl ₂	E-3200	170.00
TLA DNA Polymerase, 1,000 U, 10 mM dNTPs, 10 x Reaction Buffer with MgCl ₂	E-3201	575.00
CycleScript Reverse Transcriptase, 10,000 U, 5 x Reaction Buffer, 100 mM DTT, 10 mM dNTPs	E-3131	165.00
CycleScript Reverse Transcriptase, 50,000 U, 5 x Reaction Buffer, 100 mM DTT, 10 mM dNTPs	E-3132	540.00
5 x Reaction Buffer for CycleScript Reverse Transcriptase	E-3133	55.00
MMLV Reverse Transcriptase, 10,000 U, 5 x Reaction Buffer, 100 mM DTT	E-3121	48.00
MMLV Reverse Transcriptase, 50,000 U, 5 x Reaction Buffer, 100 mM DTT	E-3122	198.00

BIONEER, INC.

1000 Atlantic Avenue
 Alameda, CA 94501-1147 USA
 Toll Free : (877) 264 - 4300
 Fax : (510) 865 - 0350
 Email : order.usa@bioneer.us.com
 Support: support@bioneer.us.com

BIONEER CORPORATION

49-3 Munpyeong - dong, Daedeok - gu
 Daejeon 306 - 220 South Korea
 Phone: 82-42-930-8777
 Fax: 82-42-930-8688

www.bioneer.com

BIONEER
 Innovation • Value • Discovery