

# Rocketscript™ Reverse Transcriptase

● **Introduction :**

RocketScript Reverse Transcriptase is a new M-MLV originated Reverse Transcriptase optimized to synthesize first-strand cDNA from purified Poly(A) or total RNA template. RNA targets from 100 bp to > 10 kb can be detected with the RocketScript RT.

The amount of starting material can vary from 1 pg to > 1 ug of total RNA. RocketScript Reverse Transcriptase is a version of M-MLV RT that has been engineered to provide increased thermal stability.

The enzyme can be used to synthesize cDNA at a temperature range of 42 - 70 °C, providing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases.

● **Principle:**

RocketScript Reverse Transcriptase is genetically engineered, thermal stable M-MLV with enhanced thermal stability and outstanding processivity. The enzyme also features increased specificity and improved efficiency allowing efficient reverse transcription of RNA molecules with complex secondary structures

● **Contents**

Component	Quantity
RocketScript Reverse Transcriptase (10,000 U)	50 µl
5 X Reaction Buffer	0.5 ml
100 mM DTT	0.3 ml
10 mM dNTP	0.2 ml
RNase Inhibitor(100ng/ µl)	50 µl

[Cat. No.]

E-3141

[Lot No.]

[Conc]

200 U/µl

Exp.Date:

● **Protocol :**

**A. First-strand cDNA synthesis [20 µl reaction volume ]**

**Step 1** (No Incubation)

Total RNA (1 ug) or RNA (5 – 100 ng) X µl  
 Oligo dT or random primer (10 – 100 pmoles) X µl  
 water (RNase- and DNase- free) Variable volume

**Step 2**

5 X RocketScript Reaction Buffer 4 µl  
 100 mM DTT 2 µl  
 10 mM dNTP (Variable volume) or 2 µl  
 RNase inhibitor (50 – 100 ng) X µl  
 RocketScript Reverse Transcriptase (200 U/µl) 1 µl  
 DEPC-D.W Variable volume  
 Total (Step1 + Step 2) 20 µl

**B. Reverse Transcription protocol**

1. Cyclic Reverse Transcription (CRT)

Step	Temperature				Time	Cycle
	dN6	dN12	DT20	Specific Primer		
Primer annealing	15°C	30°C	37°C	T <sub>m</sub> Value	10-30 sec	12
cDNA synthesis	50°C				4 min	
Melting secondary structure & cDNA synthesis	55°C				30 sec	
Heat inactivation	95°C				5 min	1


2. Fixed Temperature Reverse Transcription (FTRT)

Step	Temperature				Time
	dN6	dN12	dT20	Specific Primer	
Primer annealing	15°C	30°C	37°C	T <sub>m</sub> Value	1 min
cDNA synthesis	42~70°C				10~60 min
Heat inactivation	95°C				5 min

**Note :** Store the buffer by low aliquots at -20 °C to minimize degradation of the DTT

Note.: Or Perform cDNA synthesis reaction as follow 50 °C 60min(cDNA synthesis), 94 °C,5min(RTase inactivation). RT reaction temperature should be selected to fit the T<sub>m</sub> value of Primers

● **Experimental Data**



**Figure 1. Amplification results of *RocketScript* Reverse Transcriptase using *myc* compared with *supplier I***

Reverse transcription condition: Incubation at each temperature 45, /50, / 55, /60, /65, /70 °C for 1 hr, inactivation at 95°C for 5 min

Primer set: human *myc* 495bp set


Lane M : 1 kb DNA Ladder

Lane 1: 10 ng Human total RNA from HeLa cell

Lane 2: 1 ng Human total RNA from HeLa cell

Lane 3: 100pg Human total RNA from HeLa cell

Lane 4: 10 pg Human total RNA from HeLa cell



**Figure 2. Comparison of long kb amplification between *RocketScript* RTase and RT enzyme from *supplier I***

RT reaction conditions were performed according to each manufacturer's recommendation.

All cDNAs were amplified with AccuPower®Hotstart PCR Premix (K-5050) from Bioneer

Note *supplier I* shows inhibition in high input concentration of total RNA

Lane 1: 1 µg Human total RNA from HeLa cell

Lane 2: 100 ng Human total RNA from HeLa cell

Lane 3: 10 ng Human total RNA from HeLa cell



**Figure 3. Comparison of amplification quality between *RocketScript* RTase and other company's RTase**

Sensitivity test. Target gene Expression Level.

Lane M : 1 kb DNA Ladder

Lane 1: 100 ng Human total RNA from HeLa cell

Lane 2: 10 ng Human total RNA from HeLa cell

Lane 3: 1 ng Human total RNA from HeLa cell

Lane 4: 100 pg Human total RNA from HeLa cell

● **Ordering Information**

Cat. No.	Description
E-3141	<i>RocketScript</i> reverse Transcriptase, 10,000 U (50 Rxn)
E-3142	<i>RocketScript</i> reverse Transcriptase, 50,000 U (250 Rxn)