

Cell-Free Synthesis of a High Molecular Weight Protein Using ExiProgen™ Automated Protein Synthesis SystemYousang Cho¹, Junho Jung¹, Nayoung Lee¹, and Jiwon Han^{1*}¹Molecular Science laboratory, Bioneer Corporation, 8-11Munpyeongseoro, Daedeok-gu, Daejeon 306-220, Korea

*Email: exiprogen@bioneer.co.kr, exiprogen-support@bioneer.co.kr

Introduction

Cell-free protein synthesis (*in vitro* protein transcription and translation) is a common method to synthesize a desired protein in a rapid and efficient manner. It is accomplished by adding template DNA (containing the coding region of the protein of interest) into a single tube which contains cell extracts and other ingredients essential for protein synthesis. Because this method does not require a separate cell line selection step, it is able to yield diverse types of protein in a very short period of time when compared to *in vivo* protein expression. Hence, cell-free protein synthesis has potential advantages in the high-throughput synthesis of proteins (1-3).

ExiProgen protein synthesis system applies automation to *in vitro* protein expression and magnetic bead-based His-Tag affinity purification and yields highly-pure proteins. *In vitro* protein synthesis using ExiProgen (Bioneer, Cat. No. A-5041) has a simple workflow; template DNA preparation, loading template DNA into ExiProgen EC1 protein synthesis kit (Bioneer, Cat. No. K-7300), loading the kit components onto the deck of ExiProgen, and starting the system. ExiProgen EC1 protein synthesis kit contains optimized *E. coli* extract which has T7 RNA polymerase and ribosomes, all other required components such as amino acids and an energy source for effective and efficient *in vitro* transcription and translation. The kit also contains Ni-NTA magnetic beads for purification of expressed his-tagged proteins. ExiProgen EC1 protein synthesis kit has 16 reaction wells so that up to 16 different kinds of highly pure proteins can be obtained at the same time within 6 hours of adding template DNA which can be in the form of a plasmid or linear PCR products.

In this study, we synthesized *Bacillus megaterium* BM3 protein *in vitro* using ExiProgen protein synthesis system to demonstrate ExiProgen's capability to automatically express and purify proteins as large as 117 kDa.

Methods and results

To synthesize BM3 protein (MW 117 kDa), 10 ug of BM3 expression vector for *in vitro* expression (a kind gift from Dr. Dong-Myung Kim at Chungnam National University) was added into a reaction well of the protein expression cartridge (cartridge 2) of ExiProgen




Figure 1. Purification and expression of BM3. ExiProgen is able to synthesize proteins as large as 120 kDa. Lane M: AccuLadder™ Protein Size Marker (Broad), Lane E: Unpurified expression sample. Lane P: Purified sample

EC1 Protein Synthesis Kit. The kit contents, including protein purification cartridge (cartridge 1), *E. coli* cell extract, elution tubes, and filter tips as well as cartridge 2, were then placed in the correct position on the deck of the ExiProgen. Next, ExiProgen was run after selecting the protocol number 902 as described in the kit manual. The run was finished in less than 6 hours and 250 µl of purified protein samples were collected in the elution tubes. To check the synthesis of the target proteins, the expression samples from J section of cartridge 2 and purification samples in the elution tubes were run in the SDS-PAGE gel. Samples were prepared for SDS-PAGE gel analysis as described in the manual. The SDS-PAGE result showed that BM3 was expressed and purified with ExiProgen (Figure 1). The amount of synthesized/purified BM3 was determined to be about 60 µg per reaction.

In conclusion, we successfully expressed and purified BM3 with ExiProgen, indicating that ExiProgen is able to synthesize proteins as large as 117 kDa. Since it is automated, easy-to-use, and provides rapid protein synthesis, ExiProgen has the potential applications several research fields including identification of protein function, protein-protein interaction study and protein structure labs as well as enzyme engineering and in Biofuel research labs. ExiProgen is proving to be a valuable tool for protein scientists.

Acknowledgement

We thank Dr. Dong-Myung Kim at Department of Fine Chemical Engineering and Chemistry, Chungnam National University for generously providing the BM3 expression vector.

References

1. D-M. Kim and J. R. Swartz. 1999. Prolonging cell-free protein synthesis with a novel ATP regeneration system. *Biotechnol. Bioeng.* 66:180-189
2. A. M. Jackson, J. Boutell, N. Cooley, and M. He. 2004. Cell-free protein synthesis for proteomics. *Brief. Funct. Genomic. Proteomic.* 2:308-319
3. J-H. Ahn, H-S. Chu, T-W. Kim, I-S. Oh, C-Y. Choi, G-H. Hahn, C-G. Park, and D-M. Kim. 2005. Cell-free synthesis of recombinant proteins from PCR-amplified genes at a comparable productivity to that of plasmid-based reactions. *Biochem. Biophys. Res. Commun.* 338:1346-1352