## 100 bp Plus DNA Ladder

DNA Molecular Weight Markers

[Cat. No.] D-1035 (1 x 500 µl) D-1036 (5 x 500 µl)

[Lot No.] 0601

[Volume] 500 µl

[Concentration] 80.7 ng/µl • **Description :** 100 bp Plus DNA Ladder is specially designed for determining the size of double strand DNA from 100 to 10,200 base pairs. The DNA Ladder consists of 19 double strand DNA fragments ranging in sizes from 100 to 1,000 bp in 100 bp increments, and additional fragments of 1,200, 1,600, 2,000, 2,961, 4,000, 5,007, 5,991, 8,000, 10,200 bp. The 500, 1,000 and 2,000 bp bands are two to three times brighter for easy identification.

• Application : Size confirmation of DNA fragments.



• Quality control assay : Agarose gel analysis shows that 19 bands are clearly distinguishable. The 500 bp, 1,000 bp and 2,000 bp band must be more intense than any other band.

Note

- The DNA ladder can be applied directly onto an agarose gel.

- Do not heat before loading.

4μt on 1.0% TBE agarose gel (Length : 8.5 cm, stained with Ethidium Bromide)

- Repeated freezing and thawing should be avoided.

## Ordering Information

[Related BIONEER DNA Molecular Weight Marker Products]

Cat. No.	Products
D-1020	25/100 bp Mixed DNA Ladder, 250 $\mu\ell$
D-1030	100 bp DNA Ladder, 250 $\mu\ell$
D-1035	100 bp Plus DNA Ladder, 500 $\mu\ell$
D-1040	1 Kb DNA Ladder, 500 $\mu\ell$
D-1050	Lambda DNA / Hind III Markers, 100 $\mu g$
D-1060	Lambda DNA / <i>Eco</i> R I Markers, 100 $\mu g$
D-1070	Lambda DNA / $\mathit{Eco} R$ I + $\mathit{Hin} d$ III Markers, 100 $\mu g$
D-2510	Lambda DNA, 250 $\mu \mathrm{g}$

Note

For research use only. Not for use in diagnostic or therapeutic procedures.