# Viral RNA PrepMate ™

## I. Description

Viral RNA PrepMate<sup>™</sup> contains all of the required reagents for easy isolation of total RNA by the guanidinium salt-lysis method. Possible sources include tissue, cultured cells, leukocytes, and serum. The extracted RNA can be used for molecular cloning, RT-PCR, Northern analysis, dot blot hybridization, poly-A (+) selection, RNase protection assays, and so on.

### II. Kit Components

Viral RNA PrepMate<sup>™</sup>(Cat.No.:K-3060)

Lysis Buffer	100 ml
RNA Hydration Solution	10 ml

### III. Protocols

A.Extraction of total RNA from tissue

- 1. Break tissue into smaller units.
- 2. Add 1ml lysis buffer to 50-100mg tissue
- 3. Homogenize the sample with a glass-Teflon homogenizer and place at 4  $^\circ\!\!C$  for 10 min.
- 4. Add 750  $\mu l$  chloroform to the sample. Mix gently and place at 4  $^\circ\!C$  for 15 min to precipitate proteins.
- 5. Centrifuge at 13,000 rpm for 20 min at  $4^{\circ}$ C.
- 6. Transfer the supernatant and add to it an equal volume of isopropanol. Store at -20  $^{\circ}$ C for one hour.
  - If the sample is placed at a low temperature (below -20°C) or 2 volumes of ethanol are added, too much salt will be precipitated.
- 7. Centrifuge at 13,000 rpm for 20 min at  $4^{\circ}$ C.
- 8. Wash pellet twice with 70% DEPC-treated ethanol.
- 9. Remove ethanol completely with cautious aspiration; spinning down may be necessary to force droplets to the tube bottom.
- 10. Remove remaining traces of ethanol by incubating open tube on a heating block for 5 min at 56  $^\circ$ C.
- 11. Resuspend the pellet in RNA hydration buffer (DEPC treated DW).
- 12. Add 6 units of DNase I and incubate for 2 hr at 37  $^\circ$ C. Inactivate DNase I by heating at 94  $^\circ$ C for 10 min.
- 13. Quantify the yield of total RNA with a spectrophotometer. Store at -70  $^\circ\!\!\mathbb{C}.$ 
  - Use RNA Hydration Buffer as a blank.
  - When OD<sub>260</sub> is 1.0, the concentration of total RNA is 40ug/ml.

#### B.Extraction of total RNA from leukocytes or cultured cells

- 1. Add 1ml lysis buffer to the suspended cells and vortex it gently in an inverted position for 10 seconds.
- 2. Incubate for 10 min at  $4^{\circ}$ C.
- 3. Follow the step III-A4 to III-A13.
- C. Extraction of total RNA from serum
- 1. Add 800ul lysis buffer to 200ul serum.
- 2. Mix carefully and place at 4  $^\circ\!\!{\rm C}$  for 10 min.
- 3. Follow the step III-A4 to III-A13. In order to proceed completed lysis, it is very important to suspend cell pellet completely.