

siRNA FAQ

Synthesis and Order

Q1. How is siRNA synthesized?

The most commonly used method in siRNA synthesizers is the 'phosphite triester' method, developed by Koster. This method builds the phosphodiester backbone of the RNA molecule by using β -cyanoethyl phosphoramidites. (Nucl. Acids Res. 1984, 12, 4539 ; Tetrahedron Lett. 1983, 24,5843) The synthesis begins with a solid support (CPG) with an attached starter nucleoside. A deblocking-coupling-capping-oxidation cycle is repeated to reach the desired length of siRNA molecule. After synthesis, ammonia is used to detach the siRNA from the solid support. The siRNA then undergoes deprotection and TBDMS treatment before being purified. The end result of these processes is pure siRNA. Each synthesized ss-RNA is annealed to produce the final duplexed siRNA form.

Q2. How do I order pre-designed siRNAs?

The *AccuTarget* Genome-wide Predesigned siRNAs are designed for Human, Mouse and Rat genomes. The desired target can be searched by Symbol, Gene No, Refseq accession Number or Description. Each query will yield three (3) candidates with the highest predicted efficiency score. You may select the desired quantity at that point. Also, you may select purification type and the guaranteed nmole quantity as well.

Q3. When can I expect delivery?

Delivery depends on the *AccuTarget* product that is ordered. Custom siRNA orders can be made and delivered in 10 business days. *AccuTarget Predesigned siRNA* can be delivered within 5 days of order receipt. Library products, orders of < 10 96-well plates can be shipped within 3 days of order receipt, and shipping dates for higher quantities above that will have to be discussed on a case-by-case basis.

Q4. Up to how many bases can Bioneer synthesize?

Bioneer offers single strand RNA synthesis service up to 50mers. (not siRNA) The same price is applied up to 30mer siRNA including overhang.

Q5. What form will my order be in?

For Genome-Wide Predesigned siRNAs, Validated siRNAs, siRNA Libraries and Control siRNAs, both the sense and antisense strands are synthesized at equimolar concentrations, verified via MALDI-TOF, then annealed and delivered in duplexed, lyophilized form. You may reconstitute the siRNA

with a buffer of your choice or with ultrapure water that we provide with every order. We recommend reconstituting to 100 μ M. For Custom siRNA orders, the order is processed by the same method as above. We recommend 50 μ M reconstitution for Custom siRNA orders. When ordering Custom siRNA you must select the "Annealing Service" to receive your order in annealed form. If you choose not to use our annealing service, you can use 1X annealing buffer and follow the annealing protocol included with your Custom siRNA order.

Q6. I want to conduct an in vitro experiment. What scale should I choose? What purification should I select?

With a 10 nmole scale siRNA order, you can transfect one hundred (100) 96-well plates at 100 nM per transfection. Unless you are planning to conduct an in vivo experiment, the Bio-RP Purification will yield outstanding results. We recommend HPLC purification for in vivo experimental use.

Q7. Can Bioneer synthesize chimeric RNA?

Yes we can. We have the ability to synthesize chimeric RNA molecules that incorporate dA, dC, dG, and dT DNA bases. We can also incorporate 2'-O-Methyl and 2'-F(rU, rC) into the RNA.

Q8. I would like to order over four (4) siRNA candidates for a single gene. How do I order?

Our Genome-wide predesigned siRNAs provide three (3) candidates per target gene. In order to request more than four (4) candidates, order via Custom siRNA request. The siRNA sequences can be verified only after the order has been submitted. If you would like to compare the sequences with publications or to modify your order, please email us

Q9. The Genome-wide predesigned siRNA didn't work like I expected. What do I do?

When purchasing 3 siRNAs for the same gene, Bioneer guarantees at least an 80% reduction in the target mRNA level for 2 of the 3 siRNAs. If there is not a > 80% reduction in the mRNA level of the target gene, Bioneer will supply 2 siRNAs free of charge* when the customer provides the experiment information.

Q10. Are phosphate groups present on the 5' or 3' ends of the synthesized siRNA?

Unless explicitly stated, the 5' and 3' ends are capped with -OH groups. Therefore, to order 5' phosphate-capped siRNAs, you must request for 5' phosphorylation modification.

Modified siRNA

Q1. What types of modified siRNAs are there?

In order to maximize siRNA applications, modified* siRNAs are available. Although there are many ways to modify siRNAs, the most convenient method is to use phosphoramidites to label siRNAs during synthesis. For a 5' siRNA modification, synthesis is first conducted as normal. At the final stage, however, a labeling phosphoramidite is attached to the 5' end, resulting in a 5' labeled siRNA. 3' labeling is somewhat complicated in that we must start with a labeled phosphoramidite attached to a solid support and build the molecule from there. You can incorporate the label of your choice to the siRNA sequence by using a labeled phosphoramidite at the 5' end of deoxyuridine, therefore allowing for an internal modified siRNA.

The following modifications are offered for our siRNAs.

- 5' Amine - 3' amine modification
- 5' Phosphorylation & 3' phosphorylation
- 5' Thiol - 3' Thiol Modification
- 5' Biotin - Internal Biotin-dT Modification
- 5' Fluorescein - 3' Fluorescein modification
- 3' TAMRA modification

*Certain items may not be available in all countries.

Handling and Storage

Q1. How to I store my siRNAs and how long do they keep?

siRNAs can normally be kept stable at -20°C for over 1 year. The lyophilized form is especially stable and has a long shelf-life. Although solution-dissolved siRNAs can be stable, contamination of the reconstitution solution with RNase will degrade the product. Also, repeated freeze-thaw cycles accelerate the degradation process. Therefore, we recommend that after you receive the siRNA stock, you reconstitute and make several aliquots to avoid such freeze-thawing. Because the phosphodiester bonds of the RNA can break under high pH conditions, we ask you to take caution, and recommend reconstituting in ultrapure water provided.

Q2. How do I reconstitute my siRNA?

If you would like to keep your siRNA in solution, we recommend reconstituting with DEPC-treated water that we provide with your order to maximize stability.

Q3. How to I store my fluorescent dye modified siRNA?

Photobleaching may occur if the fluorescent dye modified siRNA is exposed to light for prolonged periods of time. Therefore we recommend that you store such siRNAs in a dark container and store that container in a dark place.

Purification

Q1. I would like to know how the synthesized siRNA is purified.

To purify newly synthesized siRNA, we have several methods including Bio-RP and HPLC. We select a purification method depending on the end-use of the product:

1. Bio-RP: Phosphoramidites used for siRNA synthesized with the Bio-RP synthesizer have a DMT group protecting the 5' end to prevent accessory reactions from occurring during synthesis. If the DMT group is not removed at the end of the synthesis cycle ('trityl-on' mode), only the desired length N-mer siRNA will have the 5'-DMT group, while the shorter N-1, N-2 etc. synthesis by-products will not have a DMT group attached. Because the DMT group is extremely hydrophobic, it will interface very well with the RP resin. By taking advantage of this fact, and the fact that DMT groups are not present in the N-1, N-2 etc. by-products, we are able to efficiently purify the desired N-mer siRNA. Unless intended for in vivo experiments, the Bio-RP Purification method will yield excellent experimental results.

2. HPLC: For experiments requiring extremely pure siRNAs such as in vivo experiments, Bio-RP purification may be insufficiently pure. In this case, HPLC is commonly used to reach the desired purity level. We use a reversed-phase resin for purification, and this will routinely yield up to 98% purity for 21-mer siRNA molecules. We recommend selecting HPLC purification if you are intending to use the siRNA for in vivo experiment use.

Quantification and Concentration

Q1. What does the O.D. value mean?

Because the amount of synthesized siRNA is extremely small, it is very difficult to quantify by weight. Therefore, we capitalize on the light-absorption properties of nucleic acids to indirectly measure the amount of product. The amount of light absorbed by a given amount of product is expressed as an O.D. (Optical Density) value, and is directly related to the amount of product in solution.

Q2. How do I calculate the amount of siRNA from the O.D. value?

After measuring the amount of 260 nm light absorbed by the siRNA, the following formula is used to calculate the actual amount of siRNA $O.D. = \epsilon C$. The ϵ is the extinction coefficient which is unique for different materials, and the C stands for the siRNA concentration. If one knows the extinction coefficient and the O.D. value of the siRNA, the concentration can be calculated by substituting those values in the formula above. The ϵ value for an siRNA molecule is the sum of the extinction coefficient of each base constituting the molecule. The extinction coefficients for each base at 260 nm is are as follows:

rG : 11.5 ml/ μ mole

rC : 7.2 ml/ μ mole

rA : 15.4 ml/ μ mole

rU : 9.9 ml/ μ mole

Therefore, the extinction coefficient for the entire siRNA molecule can be calculated by multiplying the number of bases by the corresponding extinction coefficient and adding all four products together.

Q3. If my 18-mer siRNA (3G, 4C, 5A, 6U) O.D. value was 0.7, then how much siRNA do I have?

Let's start by calculating the extinction coefficient (ϵ) of the entire siRNA molecule. $\epsilon = 11.5 \times 3 + 7.2 \times 4 + 15.4 \times 5 + 9.9 \times 6 = 199.7$ (ml/ μ mole) Therefore, using the $O.D. = \epsilon C$ formula, we can calculate the final concentration C to be $C = 0.7 / 199.7 = 0.0035$ (μ mole/ml) = 3.5 (nmole/ml). This value is expressed as the 'total nmole' value in your copy of the oligo synthesis report.

Q4. How do I calculate the molecular weight of the synthesized siRNA molecule?

Substitute the number of each base into the following formula: $M.W. = (NA \times 329.208) + (NC \times 305.183) + (NG \times 345.207) + (NU \times 306.168) - 63.98 + 2.016$
NA = Total number of rA
NC = Total number of rC
NG = Total number of rG
NU = Total number of rU

Q5. Can I know how many ng the synthesized product is?

Normally, we will fulfill an order with a guaranteed nmole amount, and the synthesis report will also report the final amount in nmoles. If you must have the ng amount to calculate for an experiment, you can convert from nmole to ng by using the formula below. We make it easy for you by giving you the molecular weight of the siRNA sequence in the report. $Molecular\ Weight\ (g) \times mole\ (nmole) = Mass\ of\ siRNA\ (ng)$

Q6. How do I reach a target concentration?

Each synthesis report gives you a 'volume for 100 pmole/ μ l' value that stands for the volume of DEPC-treated water you need to add to achieve a 100 pmole/ μ l concentration. As an example, if the 'volume for 100 pmole/ μ l' values for two siRNAs are 100 and 200, then you can reconstitute with 100 μ l and 200 μ l of DEPC-treated water, respectively, to reach a 100 pmole/ μ l concentration.

In other words, the siRNA contained in each tube is:

$$100 \mu\text{l} \times 100 \text{ pmole}/\mu\text{l} = 10,000 \text{ pmole} = 10 \text{ nmole},$$

$$200 \mu\text{l} \times 100 \text{ pmole}/\mu\text{l} = 20,000 \text{ pmole} = 20 \text{ nmole}$$

Although the final use will dictate the concentration of siRNA, we found that for average use, 100 pmole/ μ l is ideal. That is why we provide you with DEPC-treated water volume for reconstitution to 100 pmole/ μ l.

Q7. How do I convert between molarity and moles?

The standard concentration units for oligomers is given in M (mole/L), and the prefixes such as μ (micro), n (nano), p (pico) etc. describe the scope of the unit. The following prefixes are used not only for M, but for other measurement units such as length, mass etc.

$$10^{-1} = \text{deci [d]} \quad 10^1 = \text{deca [da]}$$

$$10^{-2} = \text{centi [c]} \quad 10^2 = \text{hecto [h]}$$

$$10^{-3} = \text{milli [m]} \quad 10^3 = \text{kilo [k]}$$

$$10^{-6} = \text{micro } [\mu] \quad 10^6 = \text{mega [M]}$$

$$10^{-9} = \text{nano [n]} \quad 10^9 = \text{giga [G]}$$

$$10^{-12} = \text{pico [p]} \quad 10^{12} = \text{tera [T]}$$

$$10^{-15} = \text{femto [f]} \quad 10^{15} = \text{peta [P]}$$

$$10^{-18} = \text{atto [a]} \quad 10^{18} = \text{exa [E]}$$

$$10^{-21} = \text{zepto [z]} \quad 10^{21} = \text{zetta [Z]}$$

$$10^{-24} = \text{yocto [y]} \quad 10^{24} = \text{yotta [Y]}$$

$$1 \text{ pmole}/\mu\text{l}$$

$$= 1 \times 10^{-12} \text{ mole} / 1 \times 10^{-6} \text{ L}$$

$$= 1 \times 10^{-6} \text{ mole/L}$$

$$= 1 \mu\text{mole/L}$$

$$= 1 \mu\text{M}$$

Therefore, μ M and pmole/ μ l are one and the same.

Experiments

Q1. What are some precautions for a siRNA experiment?

First, because not all siRNAs will knock-down the target gene with identical efficiency, you should try 2-3 different sequences to find the best siRNA. Second, to make sure that the knock-down affects downstream protein expression, mRNA levels should also be measured. Thirdly, verify the knock-down phenotype by using another siRNA designed for the same target gene and show that the same phenotype appears.

Q2. This is my first siRNA experiment. How do I set my experimental conditions?

One of the most important factors in a siRNA experiment is the assessment of whether the siRNA gets delivered into the cell. Bioneer offers positive controls that can easily indicate whether the siRNA is being delivered successfully.

Q3. To what cell density should I culture my cells before siRNA transfection?

For siRNA transfection, we recommend that the cell density be approximately 60-70%. (HeLa cell, 6-well plate standard: 1.5×10^5 cells/well. We also strongly encourage user optimization of these figures)

Q4. What concentration of siRNA should I use for transfection?

We recommend a starting concentration of 100 nM (100 pmol/well), but strongly advise empirical optimization for the cell line and conditions of your experiment.

Q5. How do I use 10 nmole of siRNA to transfect cells at 100 nM?

This is a source of confusion for many people. In order to transfect a single well with 100 nM siRNA where the transfection volume is 1 ml, you need 100pmole of siRNA. 2 μ l of 50 μ M (50 pmole/ μ l) stock siRNA solution in 1 ml will yield 100 pmole. If you were to order 10 nmole of a siRNA, it will be sufficient to transfect 100 wells at 100 nM (100 pmole/ml).

Q6. What transfection reagent do I use?

Because each cell line will have different transfection efficiencies for every transfection agent, we recommend that you select the transfection reagent most suitable for your cell line.

Q7. How do I verify my siRNA transfection efficiency?

You can easily verify the transfection efficiency by transfecting your cells with NC-FITC and observing the cells with a fluorescence microscope. The NC-FITC can also be used as a test reagent to optimize the transfection concentrations of both the siRNA and the transfection reagent.

Q8. Is there a way to know my siRNA's knockdown efficiency beforehand?

The predicted siRNA knockdown efficiency can be seen at the target gene search screen by clicking the siRNA Number. Validated siRNA efficiencies will be presented in blue, while predesigned siRNA efficiencies will be displayed in grey.

Q9. What controls are used for a siRNA experiment?

The AccuTarget Negative Control siRNA is a non-targeting siRNA that has low sequence homology to all known Human, Rat and Mouse sequences. Therefore, it can be used as a convenient negative control for all Human, Rat and Mouse siRNA experiments. The AccuTarget™ Positive Control siRNAs (human) shows great knockdown efficiency for the target gene. We target GAPDH as an endogenous control gene, and also have positive control siRNAs for reporter systems such as GFP and Luciferase.

Q10. How do I verify the siRNA knockdown efficiency?

The siRNA knockdown efficiency can be verified through various techniques including qPCR, Northern Blot, Western Blot etc.

Q11. What are *AccuTarget* Genome-wide Predesigned siRNAs?

These siRNAs were predesigned for your convenience using our Turbo si-Designer, a siRNA design algorithm co-developed by us and KRIBB (Korea Research Institute of Bioscience & Biotechnology). Using this revolutionary algorithm, we have constructed a high-knockdown efficiency *AccuTarget* Genome-wide Predesigned siRNA database. As of this writing, our website contains predesigned siRNAs for Human (17,842 genes), Mouse (17,117 genes), and Rat (9,392 genes).

Q12. What are *AccuTarget* validated siRNAs?

After synthesizing siRNAs using our Turbo si-Designer, a siRNA design algorithm co-developed by us and KRIBB (Korea Research Institute of Bioscience & Biotechnology), *AccuTarget* Validated siRNAs have been measured for knockdown efficiency using qPCR. We verified that each siRNA shows at least a 70% knockdown efficiency. *AccuTarget* Validated siRNAs can be ordered with their corresponding real-time qPCR primer sets. Also, if you order 10 or more siRNAs, they will be delivered in flexible siRNA Library (tube or plate type) format.

Q13. Can I order siRNAs targeted for a specific gene?

You may order siRNAs for the target gene of your choice by entering the gene name or Accession number at our website and selecting one of the listed siRNA candidates. You may also order custom designed siRNA if your model organism is not Human, Mouse or Rat.

Q14. What are precautions for handling siRNA?

- Wear gloves and a mask when handling siRNA.
- Always use RNase free tubes and pipette tips.
- Only use DEPC treated water when diluting siRNA.
- Store fluorescence-labeled siRNA in a dark location.